

# vPad-A1™

Datrend Systems' vPad-A1, based on our revolutionary Vision-Pad Technology™, is an "all-in-one" patient simulation system. vPad-A1 is modular and is comprised of a Multi-Parameter Patient Simulator, SpO2 test module, and a Non-Invasive Blood Pressure simulation module which may be used together, independently, or in various combinations. An Android handheld device or vPad tablet provides the user interface.

## Key Features:

- 12 lead ECG simulation
- ST Segments: 8 elevated and 8 depressed
- Axis Deviation: Normal (intermediate), horizontal, and vertical.
- Neonatal Mode
- ECG Performance Testing
- Over 60 Arrhythmia selections
- 2 channel Invasive BP simulation
- Temperature and Respiration simulations
- Pacer simulations
- Cardiac output
- SpO2 Pulse Oximeter Simulator and Non-Invasive BP modules, compatible with all major industry manufacturers
- Auto Settings
- Auto Sequences
- Test reports - with user entered results
- Control via Bluetooth or USB



Innovation by design

## 宇宏企業股份有限公司

台北：三重區重新路五段 609巷 16號 5樓之7

TEL:02-22782269 FAX:02-22782059

新竹：新竹縣竹北市台元 26 號 4 樓之 10

TEL:03-552-6779 FAX:03-552-6739

[http:// www.radiotek.com.tw/](http://www.radiotek.com.tw/)

E-mail : [salesmag1@radiotek.com.tw](mailto:salesmag1@radiotek.com.tw)

## vPad-A1™

The vPad-A1 contains the following modules: A1 Base module, vPad-PS, vPad-O2, and vPad-BP.

The vPad-A1 App is an application developed for an Android handheld device, which serves as the user interface for the vPad-A1 system. It is designed using the same principles as all other Datrend vPad Apps, which are intuitive, ergonomic, and customizable.

The vPad-A1 Base serves as a communication bridge between the Android tablet and the rest of vPad-A1 units via Bluetooth or USB. A DACOM bus connection provides an interface to other vPad devices. The Base also provides power to vPad-PS and vPad-O2.

vPad-PS provides six of the eight available simulations: ECG, two-channel IBP (invasive blood pressure), respiration, temperature, and cardiac output.

vPad-O2 provides SpO2 simulation for testing pulse oximeters, and is compatible with all major pulse oximeter manufacturers.

vPad-BP provides NIBP (non-invasive blood pressure) simulation supporting all major NIBP manufacturers.



130 - 4020 Viking Way, Richmond BC, Canada V6V 2L4  
Phone: 604.291.7747 • Fax: 604.294.2355  
Toll-free (North America only): 800.667.6557  
Email: customerservice@datrend.com

© DATREND SYSTEMS INC. October 2016

### ECG General:

Full 12-Lead ECG; independent outputs for each signal lead  
 - color coded to AHA and IEC Standards.  
 Output Impedances: 500, 1000, 1500, & 2000 ohms  
 ECG Amplitude: 0.05 - 5.5 mV  
 Amplitude Accuracy:  $\pm$  (2% setting + 0.05 mV)  
 High Level ECG: 500x lead II signal  
 High Level Accuracy:  $\pm$  5%  
 Rate Accuracy:  $\pm$  0.25 BPM

### Normal Sinus Rhythm:

Rates: 10-360 BPM, 1 BPM steps, Accuracy  $\pm$  0.25BPM  
 user defined presets (15), user input specific rates  
 Amplitudes (Lead II): 0.05mV to 0.5mV (0.05mV steps;  
 0.5mV to 5.5 mV (0.25mV steps)  
 Neonatal Mode: ECG QRS width is reduced from 80ms  
 to 40ms.  
 Artifact: 50Hz, 60Hz, muscle, baseline, respiration  
 Axis Deviation: Normal, horizontal, and vertical.

### ECG Performance Testing:

Square Wave: 0.125, 2, 2.5Hz  
 Triangle Wave: 0.125, 2, 2.5Hz  
 Pulse: 30, 60 BPM with 60ms pulse  
 Sine Waves: 0.05 - 200 Hz.  
 QRS and R Wave Detection Test:  
 Rate: 30 - 250 BPM triangle wave  
 Width: 8 - 200ms  
 ST Segment Adjustment (Lead II):  
 Rate: 60 BPM; ST Segment:  $\pm$  80% of ECG amplitude  
 Tall T wave:  
 Rate: 80 BPM; ST Segment: 0 - 150% of ECG amplitude

### Fetal / IUP(ch1 only) Simulations:

Fetal heart rates: 60 to 240 BPM 1 BPM steps  
 12 Preset rates, user defineable  
 Uniform, Early and Late Deceleration,  
 Uniform Acceleration  
 Dynamic intrauterine pressure (IUP) waveform:  
 Positive bell shaped pressure curve  
 Peak pressure: 50 or 90 mmHg,  
 Contraction duration: 90 sec  
 IUP Period: 2, 3, 5 min and Manual  
 Pressure transducer sensitivity: 5 or 40 m v/v/mmHg  
 Input/output impedance: 300 ohms  $\pm$ 10%

### 2 Blood Pressure Channels:

Electrically Isolated Channels  
 Transducer Sensitivity: 5 or 40  $\mu$ V/V/mmHg  
 Input/output impedance: 300 ohms  $\pm$ 10%  
 Excitation : 2 to 16 Vp; DC to 5000Hz  
 Calibrated Rate: 80 BPM normal sinus rhythm  
 Static Levels BP1/2:  
 -10 to 400 mmHg in 1 mmHg steps  
 15 User defined presets; user input specific pressures  
 Accuracy:  $\pm$  (1% of setting + 1mmHg)

### Dynamic Simulations:

Arterial (120/80)  
 Arterial (90/40)  
 Arterial (160/110)  
 Radial Artery (120/80)  
 Left Ventricle (120/0)  
 Right Ventricle (25/0)  
 Pulmonary Artery (25/10)  
 Pulmonary Artery Wedge(25/2)  
 Right Atrium [CVP] (120/0)  
 Left Atrium (14/4)  
 Swan-Ganz (channel 1 only)  
 Automatic (every 15, 25sec) with Pause  
 Manual, advance is manually triggered  
 Artifact/Respiration (larger of):  
 5mmHg or 5%  
 10mmHg or 10%

### Pacemaker:

Pulse Amplitude: -700mV to +700mV  
 Pulse Polarity: Positive or negative.  
 Pulse Width: 0.1, 0.2, 0.5, 1.0, 2.0 ms  
 Accuracy :  $\pm$ (5% setting + 0.2mV) Lead II  
 Pacer Rhythm:  
 Ventricular  
 Asynchronous 75 BPM  
 Demand with frequent sinus beat  
 Demand with occasional sinus beat  
 A-V sequential  
 Non-capture  
 Non-function  
 Atrial  
 Atrial 80 BPM  
 A-V sequential

### Temperature:

20 - 42°C in 0.5°C increments  
 Accuracy:  $\pm$ 0.01 °C high precision simulations  
 (30, 32, 35, 37, 40, 42 °C)  
 $\pm$ 0.03 °C general  
 Probe Compatibility: 400 or 700 series YSI

### Respiration:

Baseline Impedance:  
 500, 1000, 1500, 2000 ohms on LEADS I, II, III  
 Impedance Variations (Delta):  
 0.05 to 1.0 $\Omega$  in 0.05 $\Omega$  increments;  
 1.0 to 5.0 $\Omega$  in 0.25 $\Omega$  increments;  
 Rates: 10 to 150 BrPM; 1 BrPM steps; 0 BrPM for APNEA  
 Apnea Selections: 12, 22, 32 seconds, and continuous  
 Respiratory Effort (Inspiration/Expiration Ratio:) 1/1, 1/2,  
 1/3, 1/4, 1/5  
 Ventilated 1/1  
 Respiration Lead LA or LL

### Cardiac Output:

Baseline Temperature: 36, 37 and 38°C,  $\pm 0.03$  °C  
 8 Inject Temperatures 0, 2, 20 & 24°C; Spacelabs and  
 Phillips  
 1 user adjustable  
 Simulations:  
 C.O. of 3, 4, 5, 6, 7l/min  
 Slow Injectate Curve  
 Faulty Injectate Curve  
 Left to Right Shunt Curve  
 Cal Pulse: 1°C for 1 second

### Arrhythmia Selections:

#### General 1

Asystole 1  
 Asystole 2  
 Asystole 3  
 PVC1 Bigeminy  
 PVC1 Trigeminy  
 PVC2 Bigeminy  
 PVC2 Trigeminy  
 Premature Atrial Contraction (PAC)  
 Nodal Premature Nodal Contraction (PNC)  
 Multifocal PVC (once)  
 Frequent Multifocal PVCs

#### Ventricular Arrhythmia (PVC1\left or 2\right)

PVC Ventricular (once)  
 PVC Ventricular (every 10th beat)  
 PVC Early, Ventricular  
 PVC R-on-T, Ventricular  
 PVC 6/Minute  
 PVC 12/Minute  
 PVC 24/Minute  
 Pair PVCs (1 time event)  
 Run 5 PVCs (1 time event)  
 Run 11 PVCs (1 time event)

### Conduction Defects:

First Degree Heart Block  
 Mobitz I, Second Degree Heart Block  
 Mobitz II, Second Degree Heart Block  
 Third Degree Heart Block  
 Right Bundle Branch Block  
 Left Bundle Branch Block

#### Fibrillations

Coarse Atrial Fibrillation  
 Fine Atrial Fibrillation  
 Coarse Ventricular Fibrillation  
 Fine Ventricular Fibrillation

### Supraventricular Arrhythmia

Atrial Tachycardia  
 Paroxysmal Atrial Tachycardia  
 Supraventricular Rhythm @ 90 & 120 BPM  
 Supraventricular Tachycardia @ 140, 150, 160,  
 180, 190, 200, 210 & 220 BPM  
 NSR @ 160 BPM

#### General 2

Atrial Flutter  
 Sinus Arrhythmia  
 Missed Beat @ 80 BPM (1 time event)  
 Miss every 10th @ 80 BPM  
 Miss every 10th @ 120 BPM  
 Nodal Rhythm  
 Sinus Bradycardia <60 BPM

#### AutoSettings

Unlimited number of user programmable,  
 simulation parameter setups available.

### Communication / User Interface:

via vPad-A1 Base Unit  
 Android 5" tablet:  
 Touchscreen User Interface  
 Wired (USB) or Bluetooth mode  
 WiFi  
 16 GB memory  
 Dual XBUS for Datrend test automation

### Power Supply:

via vPad-A1 Base Unit  
 External AC adapter  
 Internal rechargeable Li-Ion batteries (for 10 hrs  
 of simulation with full charge)

### Dimensions:

98mm x 208mm x 56mm (3.85" x 8.2" x 2.21")  
 PS Unit (incl. A1 Base)

### Weight:

660g (1.44lb) PS Unit (incl. A1 Base)  
 200g (0.44lb) wireless tablet interface

### Environment:

15°C to 40°C, 10% to 90% RH, Indoor Use Only,  
 Category II

*All specifications subject to change without notice.*



130 - 4020 Viking Way, Richmond BC, Canada V6V 2L4  
 Phone: 604.291.7747 • Fax: 604.294.2355  
 Toll-free (North America only): 800.667.6557  
 Email: customerservice@datrend.com © DATREND SYSTEMS INC. SEPTEMBER 2016

**Saturation (SpO<sub>2</sub>):**

Range: 30% to 100%

Increments: 1%

Presets: 6, user definable

Range of adjustment and presets may vary according to pulse oximeter specifications

**SpO<sub>2</sub> Accuracy:**

Saturation within DUT specified range:

±1 count + specified accuracy of the DUT

**Heart Rate:**

Range: 20 to 300 BPM

Increments: 1 BPM

Presets : 6, user definable

Accuracy: ± 0.25 BPM (sync mode)

otherwise, ±1 BPM

**Pulse Amplitude:**

Range: 0 to 100%

Increments: 1% steps.

Presets : 6, user definable

Accuracy: ± 1%

**Signal Artifact:**

Four preset simulations:

Movement

Tapping (Spike artifact)

Shivering (Tremor artifact)

Shake Table (2.5Hz Sinewave)

**Auto Presets:**

Unlimited preset patient simulations

Default Auto Presets:

Normal Adult

Hypoxia

Movement Artifact

Tachycardia

Bardycardia

Neonate

Low Perfusion

No perfusion

Tremor (Shivering Artifact)

**Alarm Tests:**

Automated test sequences for determining oximeter alarm response time to:

Low Saturation

Low Heart Rate

High Heart Rate

Low Perfusion

Signal Artifact

Five defaults, plus unlimited programmable alarm sequences

**Communication / User Interface:**

via vPad-A1 Base Unit

Android 5" tablet:

Touchscreen User Interface

Wired (USB) or Bluetooth mode

WiFi

16 GB memory

Dual XBUS for Datrend test automation

**Power Supply:**

via vPad-A1 Base Unit

External AC adapter

Internal rechargeable Li-Ion batteries (for 10 hrs of simulation with full charge)

**Dimensions:**

98mm x 208mm x 30mm (3.85" x 8.2" x 1.18") A1 Base

90mm x 160mm x 24mm (3.54" x 6.3" x 0.95") SpO<sub>2</sub> Probe

**Weight:**

440g (0.96lb) A1 base

122g (0.27lb) SpO<sub>2</sub> Probe

**Environment:**

15°C to 40°C, 10% to 90% RH, Indoor Use Only, Category II

*All specifications subject to change without notice.*

**Manufacturer Envelopes:**

GE Dinamap and Dinamap Pro  
Critikon Dinamap Plus  
Datascope Passport  
Welch-Allyn Spot LXi and 52000  
Fukuda Dynascope  
Colin Press-Mate  
...and more custom simulations

**Pressure Units:**

mmHg, mbar, kPa, inH<sub>2</sub>O, cmH<sub>2</sub>O

**BP Presets:**

Systolic/Diastolic (mmHg)  
Adult Neonatal  
255/195 150/100  
200/150 120/80  
150/100 100/65  
120/80 80/50  
100/65 60/30  
80/50 35/15  
60/30

**BP Simulation:**

Simulation type: oscillometric  
Rate Range: 20 – 240 BPM  
Accuracy: ±0.25 BPM in sync mode  
±1 BPM otherwise  
Amplitude: 0 – 2 mL  
1.25 mmHg into 500ml cuff  
Amplitude Range: 0 – 150%  
Amplitude Accuracy: better than 0.5%  
AutoSettings: unlimited, user definable

**Envelope Shift:**

± 50 mmHg max  
Minimum Diastolic: 15 mmHg  
Maximum Systolic: 275 mmHg

**Manometer:**

Pressure Range: 0.0 to 400.0 mmHg  
Accuracy: ± 0.5 mmHg  
Resolution: 0.1 mmHg

**Regulated Pressure Source:**

Pressure Range: 10.0 to 400.0 mmHg  
Accuracy: ± 0.5 mmHg  
Resolution: 0.1 mmHg

**Leak Test:**

Automatic/manual Inflation  
Automatic Timer  
Leak Test Time: 30 - 600 seconds (user defined)  
Target Pressure: 20 to 400 mmHg  
Range: 0 to 200 mmHg/min  
User Definable Presets: 12  
User Definable AutoSettings: unlimited

**OverPressure Test:**

Automatic/manual Inflation  
Range: 20 - 400 mmHg  
Release Time: 1-999 sec  
User Definable Presets: 12  
User Definable AutoSettings: unlimited

**Standard Features / Accessories:**

- Autosequences
- Unlimited User Defined Settings
- vPad-A1 Power Base / Display
- Universal Hose Adapter Kit

**Communication / User Interface:**

via vPad-A1 Base Unit  
Android 5" tablet:  
Touchscreen User Interface  
Wired (USB) or Bluetooth mode  
WiFi  
16 GB memory  
Dual XBUS for Datrend test automation

**Power Supply:**

External AC adapter  
Internal rechargeable Li-Ion batteries (for 200+ simulations with full charge)

**Dimensions:**

98mm x 275mm x 97mm (3.87" x 10.82" x 3.80")  
BP Unit (incl. A1 Base)

**Weight:**

1080g (2.38lb) BP Unit (incl. A1 Base)  
200g (0.44lb) wireless tablet interface

**Environment:**

15°C to 40°C, 10% to 90% RH,  
Indoor Use Only, Category II

*Please contact the factory for the availability of other calibration tables, or visit our web site for updates at [www.datrend.com](http://www.datrend.com)*

*All specifications subject to change without notice.*



130 - 4020 Viking Way, Richmond BC, Canada V6V 2L4  
Phone: 604.291.7747 • Fax: 604.294.2355  
Toll-free (North America only): 800.667.6557  
Email: [customerservice@datrend.com](mailto:customerservice@datrend.com) © DATREND SYSTEMS INC. SEPTEMBER 2016



**宇宏企業股份有限公司**

台北：三重區重新路五段 609巷 16號 5樓之7  
TEL:02-22782269 FAX:02-22782059  
新竹：新竹縣竹北市台元 26 號 4 樓之 10  
TEL:03-552-6779 FAX:03-552-6739  
[http:// www.radiotek.com.tw/](http://www.radiotek.com.tw/)  
E-mail : [salesmag1@radiotek.com.tw](mailto:salesmag1@radiotek.com.tw)